



上海科技大学
ShanghaiTech University

Lecture 9: Introduction to CT

Dr. Guohua Cao (曹国华)

caogh@shanghaitech.edu.cn

School of Biomedical Engineering

ShanghaiTech University

4/2/2024

Lecture 9: Introduction to CT

- What is CT
- Local contrast in CT images
- Line integral in CT
- Image pre-processing
- Sinogram
- CT Number and Hounsfield unit
- First CT
- CT Generations
- CT Hardware
 - Detectors
 - Gantry
 - Slip Ring
- Spiral CT
- Dual Source CT/Dual Energy CT
- Spectral CT


Lecture 9: Introduction to CT

- What is CT
- Local contrast in CT images
- Line integral in CT
- Image pre-processing
- Sinogram
- CT Number and Hounsfield unit

- First CT
- CT Generations
- CT Hardware
 - Detectors
 - Gantry
 - Slip Ring
- Spiral CT
- Dual Source CT/Dual Energy CT
- Spectral CT


What is X-ray Computed Tomography (CT)?

1. Produce X-rays with X-ray tube
2. X-rays pass through patient
3. X-ray detected on the other side
4. Repeat from all angles surrounding patient
5. Reconstruct cross sectional images
6. Each voxel values represent x-ray attenuation of the tissues at this volume element.



CT example images


2D slice images



Volume rendered image




Physics of CT




Principles of CT


What signals do we measure?



$$I = I_0 e^{-\sum \mu_i \Delta x} \rightarrow S = 1/\Delta x \log I_0/I = \sum \mu_i$$



How do we reconstruct CT images?



Solve a system of linear equations,

$$S_1, S_2, S_3, S_4, \dots$$


to find unknowns,

$$\mu_1, \mu_2, \mu_3, \mu_4, \dots$$

Contrast in X-ray projection image

Attenuation signal in radiography represents a summation of individual attenuation by various tissues along the x-ray beam path from x-ray focus to detector. That is, each pixel in x-ray radiographs indicates the integral of the attenuation along a line (i.e. “line integral”).

Contrast in x-ray radiographs is dominated by tissues with high attenuation such as bone and contrast media, or by thick objects. Those soft tissues with low attenuation are therefore mostly hidden.




“Local” contrast in CT slice image

In each CT slice image, contrast is given by the **difference in attenuation values of neighboring voxels**.


Contrast is determined **locally** by the tissue composition within the volume elements. That is, contrast in CT images is “local”.

This is fundamentally different from radiography!



This lung tumor (arrow indicated) would be difficult to spot from conventional x-ray radiography.

Raw signals in CT imaging



$$I = I_0 e^{-\sum \mu_i \Delta x} \rightarrow S = \log I_0 / I = \Delta x \sum \mu_i$$

I & I_0 are the raw signals/data measured by the detector.


S is the signal pre-processed from the raw data.

S is total x-ray attenuation along the beam path (which is a pencil-like line – pencil beam).


➔ CT is an x-ray attenuation based imaging modality (just like X-ray radiography).

The CT signal is called the **line integral**!



Monochromatic radiation passes through homogeneous object




Monochromatic radiation passes through inhomogeneous object



Data Pre-Processing Example




Projection as forward process – Forward Projection




Projection – a forward process,
i.e. from object to projections.

Reconstruction as backward process – back-projection is reconstruction




Reconstruction – a backward process that convert the projection images taken at different angles into a volumetric stacks of images (i.e. the object).

Reconstruction: Projections \rightarrow Sinograms \rightarrow Images




A sinogram example

- Object: a (2D) phantom containing a square, a circle, a triangle in a dark background



The sinogram is a **composite** representation of projection data of the entire 2D object.

Sinogram generation using Matlab



Matlab Code for previous slide

```
figure(33);
ph = zeros(128);
ph(64,64)=1;
ph(32,32)=1;
subplot(2,2,1);
imagesc(ph);

theta = 0:180;
[P,xp] = radon(ph,theta);
subplot(2,2,2);
imagesc(P);
axis normal;
title('Parallel-Beam Projections of 2 dots');
xlabel('\theta (degrees)');
ylabel('t');
colormap(hot);
colorbar;



ph = phantom(128); /* try P = phantom('Modified Shepp-Logan',256); for practice */
subplot(2,2,3);
imagesc(ph);
subplot(2,2,4);
[P,xp] = radon(ph,theta);
imagesc(P);
axis normal;
title('Parallel-Beam Projections of Shepp-Logan Phantom');
xlabel('\theta (degrees)');
ylabel('t');
colormap(hot);
colorbar;
```

CT Number & Hounsfield Unit

After reconstruction, we got a 'map' of attenuation coefficients at each voxel.

The attenuation coefficient at each voxel in CT images, when expressed in a special unit (i.e. **Hounsfield unit**), is called CT number.

In Hounsfield unit (HU), CT number for water voxel is 0 HU and the attenuation for air voxel is -1000 HU. CT number for all other type of material voxels is calculated as




Lecture 9: Introduction to CT


- What is CT
- Local contrast in CT images
- Line integral in CT
- Image pre-processing
- Sinogram
- CT Number and Hounsfield unit

- First CT
- CT Generations
- CT Hardware
 - Detectors
 - Gantry
 - Slip Ring
- Spiral CT
- Dual Source CT/Dual Energy CT
- Spectral CT

First Generation CT Scanner



pencil beam (1970)




(a) 1st generation: translation/rotation



Sir Godfrey Hounsfield


Typical guy engineer!

How does the first CT scanner work?



First Generation – Pencil Beam

pencil beam (1970)




(a) 1st generation: translation/rotation

Second Generation – First clinical CT scanner


Second Generation –

Partial Fan Beam

partial fan beam (1972)




(b) 2nd generation: translation/rotation



*The first clinical CT scan: Atkinson
Morley's Hospital, October 1971
England*

Third-generation CT scanner


Fan Beam
fan beam (1976)



(c) 3rd generation: continuous rotation




1980




2000




“Slice race” in Multi-Detector row CT (MDCT)




Toshiba Aquilion 320, 320 detector rows, 0.275s rotation, 0.5 mm slice thickness, introduced 2010s




Multi-slice spiral CT




Five Generations of CT




1st Generation




2nd Generation




3rd Generation



4th Generation




5th Generation




Spiral CT

CT detectors




GE's proprietary Gemstone Detector




Siemens's
proprietary
ultra fast
ceramics
detector

Single row CT detector



64-slice detector designs for multi-slice spiral CT



Diagrams of various 64-slice detector designs (in z-direction). Lee W. Goldman J. Nucl. Med. Technol. 2008;36:57-68

Rotating gantry





Gantry revolution time	Centrifugal force on tube
1s	3g
0.3s	30g

*space shutter liftoff ~ 3g



>3 tons, 0.275s per revolution!

Slip ring



A Look Inside a Slip Ring CT

X Ray
Tube

Detector
Array

Slip Ring



Note:
how most
of the
electronics
is
placed on
the rotating
gantry



18: Optimization of Protection in CT Scanner

13

Advances in CT

Table 1. Performance characteristics^a of CT in a comparison from 1972 to 2005. Significant improvements are given due to technological advances for almost all parameters. Contrast resolution reached a stable level since detector quantum efficiency was close to 100% early on.

	1972	1980	1990	2004	2005 (DSCT)
Rotation time (s)	300	5–10	1–2	0.33–0.5	0.33
Data per 360° scan (MB)	0.058	1	1–2	10–100	20–200
Data per spiral scan (MB)	–	–	24–48	200–4000	200–8000
Image matrix ^b	80 × 80	256 × 256	512 × 512	512 × 512	512 × 512
Power (kW)	2	10	40	60–100	2 × 80
Slice thickness, mm	13	2–10	1–10	0.5–1	0.5–1
Spatial resolution (LP cm ⁻¹)	3	8–12	10–15	12–25	12–25
Contrast resolution	5 mm/5 HU/ 50 mGy	3 mm/3 HU/ 30 mGy	3 mm/3 HU/ 30 mGy	3 mm/3 HU/ 30 mGy	3 mm/3 HU/ 30 mGy

^a Typical values for high performance scanners.

^b Values refer to the calculated matrix. Monitor displays often use 1024 × 1024 matrices by means of interpolation.

Kalender, W. "X-ray computed tomography", Physics in medicine and biology 51(13) (2006).


Dual source / dual energy CT



Dual energies:

140 keV and 80 keV

Decomposed into iodine, bone and tissue




Dual-energy x-rays / CT




Photon-counting CT



Conventional CT




Photon-counting CT



Photon-counting detector for color CT



Spectral CT



J. P. Schlomka, et al, "Experimental feasibility of multi-energy photon-counting K-edge imaging in pre-clinical computed tomography," *Phys. Med. Biol.*, vol. 53, no. 15, pp. 4031–47, Aug. 2008.